

Illumina Sequencing Overview: Library Prep to Data Analysis

Noel Lenny, PhD
Sr. Clinical Field Applications Scientist
28-SEP-2020



QB7845

For Research Use Only. Not for use in diagnostic procedures.

© 2019 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiniSeq, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, veri, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the US and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

illumina


Session Objectives

- **By the end of this training, you will be able to:**
 - List the major steps in the Illumina sequencing workflow
 - Library Preparation
 - Cluster Generation
 - Sequencing
 - Data Analysis
 - Discuss the sequencing by synthesis process
 - 4-Channel Chemistry
 - 2-Channel Chemistry
 - 1-Channel Chemistry

Illumina Sequencing Workflow

1

Library Preparation





3

Sequencing

STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
ST AGGCTAGGTTTCATGCTA
STA GCTAGGTTTCATGCTA
STAA CTAGGTTTCATGCTA

Library Prep is Critical for Successful Sequencing




For clustering:
Libraries must have P5 and P7 binding regions on either end of a library

For sequencing:
Libraries must have sequencing primer binding regions

For mixing samples:
Libraries must have a unique index or barcodes sequence

Illumina Sequencing Workflow

Cluster Generation




Sequencing

STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
ST AAGGCTAGGTTTCATGCTA
ST AGGCTAGGTTTCATGCTA
STA GCTAGGTTTCATGCTA
STAA CTAGGTTTCATGCTA

What is a Cluster?



Clusters are a group of DNA strands positioned closely together

Each cluster represents thousands of copies of the same DNA strand in a 1–2 micron spot



An image of fluorescently labelled clusters on a flow cell

**Single
DNA
Library**





**Amplified
Clonal
Cluster**

What is a flow cell?

Cluster generation occurs on a flow cell

A flow cell is a thick glass slide with channels or lanes

Each lane is coated with a lawn of oligos complementary to library adapters




Flow cell Architecture

Random vs Patterned


Random Flow Cell

- HiSeq™ 2500, MiSeq™, NextSeq™, MiniSeq™
- Randomly spaced clusters
- Variable Insert Sizes
- Lower Duplication Rates




Patterned Flow Cell

- HiSeq 3K/4K/X, NovaSeq™ 6000, iSeq™ 100
- Defined size and spacing
- Increased Cluster density
- Simplified imaging




Patterned flow cells


Complete control of pitch & feature size






Pitch (Center to Center)



Decrease Pitch



Decrease Feature Size




Random

Patterned

Rigid registration reduces time by skipping template generation

Traditional Cluster Generation



Hybridize Fragment & Extend


Single-stranded DNA libraries are hybridized to primer lawn

Bound libraries are then extended by polymerases

Surface of flow cell coated with a lawn of oligo pairs

Single-stranded DNA library template molecule

3' extension




Denature Double-Stranded DNA


Double-stranded molecule is denatured

Original template washed away


Newly synthesized strand is covalently attached to flow cell surface



Single-Stranded DNA




NOTE:
Single molecules
bind to flow cell in a
random pattern



Bridge Amplification


Single-stranded molecule flips over and forms a bridge by hybridizing to adjacent, complementary primer

Hybridized primer extends by polymerases



Bridge Amplification

Double-stranded bridge is formed




Denature Double-Stranded Bridge

Double-stranded bridge is
denatured

Result:


Two copies of covalently bound
single-stranded templates



Bridge Amplification


Single-stranded molecules flip over to hybridize to adjacent primers

Hybridized primer extends by polymerase




Bridge Amplification

Bridge amplification cycle is repeated until multiple bridges are formed




Linearization

dsDNA bridges are
denatured



Reverse Strand Cleavage


Reverse strands are cleaved and washed away, leaving a cluster with forward strands only



Read 1 Primer Hybridization


Sequencing primer is hybridized to Read 1 sequencing primer binding site

Sequencing primer




Patterned Flow Cell and ExAmp Technology

- ExAmp technology creates clonal clusters in each well from individual library molecules




The background of the slide is a repeating pattern of small, light-gray circular microwells. Each microwell contains a cluster of small, vertical, colored lines representing oligonucleotides (oligos) in yellow, green, and blue. These are arranged in a grid-like fashion across the entire slide.


Oligos present inside patterned
microwells




Oligos are not present on surface
between wells to control cluster size




Single stranded DNA is introduced to the flow cell, which hybridizes to an oligo present in the microwell




When hybridization occurs amplification immediately begins




Through a proprietary process, double stranded ends of the DNA are denatured and hybridize to the flow cell




Amplification occurs on both strands simultaneously




Amplification displaces double stranded
DNA




Amplification continues until fully double stranded DNA is formed




After double stranded DNA is formed, the process repeats itself




This process continues until no unused oligos are left on the flow cell surface



To prepare the strands for sequencing,
linearization is performed




Reverse strands are cleaved



Sequencing primer is hybridized to Read 1 sequencing primer binding site

Illumina Sequencing Workflow

Sequencing



STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
STAAGGCTAGGTTTCATGCTA
ST AAGGCTAGGTTTCATGCTA
ST AGGCTAGGTTTCATGCTA
STA GCTAGGTTTCATGCTA
STAA CTAGGTTTCATGCTA

Illumina Sequencing Systems

Focused Power



iSeq™ 100



MiniSeq™



MiSeq™



NextSeq™ 500/550

Production Power



HiSeq™ 2500



HiSeq™ 3000/4000



HiSeq™ X Five/Ten




NovaSeq™ 6000

Population Power

4-Channel SBS Chemistry:


HiSeq, MiSeq



Each of the four DNA bases emits an intensity of a unique wavelength


Collects four images:

- During each cycle, each cluster appears in only one of four images



A Closer Look At 4-Dye Chemistry



4-channel chemistry



2-Channel SBS Chemistry:

NextSeq 550, MiniSeq, NovaSeq 6000


- 2-channel SBS uses two images:
 - Clusters appearing in green only are **T**
 - Clusters appearing in red only are **C**
 - Clusters appearing in both images are **A**
 - Clusters not present/dark are **G**
-
- After imaging, cluster intensities are plotted and bases called accordingly



[Illumina Two-Channel SBS Sequencing Technology Technote](#)

A Closer Look At 2-Dye Chemistry

2-channel chemistry



Add 3 FI-NTP's &
nonFI-NTP + Polymerase

Incorporated FI-NTPs excited
then imaged (2 images/cycle)

Terminator & fluorescent
dye cleaved from FI-NTP

X # of bp in read

1-Channel SBS Chemistry:

iSeq 100

SBS chemistry combined with Complementary Metal-Oxide-Semiconductor (CMOS) technology

- The system uses a patterned flow cell with nanowells fabricated over a CMOS chip
- Each sequencing cycle has two chemistry steps in order to determine bases
- Two images are captured within one cycle of sequencing run


Based on the signal pattern across two images, base calls can be determined

- Intensities extracted from one image and compared to a second image result in four distinct populations, each corresponding to a nucleotide.


A Closer Look At 1-Dye Chemistry

1-channel chemistry


One Sequencing Cycle




Incorporation



Imaging



Chemistry



Imaging

Incorporation

First Image
captured

Second Chemistry Step


Second Image
captured

X # of bp in read

Sequencing by synthesis with CMOS detection









1-channel chemistry








Base	Image 1	Image 2
T	ON	ON
A	ON	OFF
C	OFF	ON
G	OFF	OFF










- The iSeq 100 System uses 1-dye sequencing, which requires one dye and two images to encode data for the four bases

Illumina Chemistry Comparison

4-Channel Chemistry				
	 A	 G	 T	 C
Image 1				
Image 2				
Image 3				
Image 4				
Result	A	G	T	C

2-Channel Chemistry				
	 A	G	 T	 C
Image 1				
Image 2				
Result	A	G	T	C

1-Channel Chemistry				
	 A	G	 T	 C
Image 1				
Image 2				
Result	A	G	T	C

----- Intermediate chemistry step

4-channel SBS

- Bases are identified using four different fluorescent dyes, one for each base and four images per sequencing cycle

2-channel SBS

- Simplified nucleotide detection by using two fluorescent dyes and two images to determine all four base calls

1-channel SBS


- Base calling uses one fluorescent dye and two images, with chemistry step in between, to determine all four base calls

Paired-End Sequencing



Paired-End Sequencing


The read product is stripped off




Paired-End Sequencing

Single-stranded template loops over to form a bridge by hybridizing with a lawn primer

3'-ends of lawn primer are extended




Paired-End Sequencing




Paired-End Sequencing

Bridges are linearized and the original forward template is cleaved




Paired-End Sequencing

Sequencing primer is hybridized to the Read 2 sequencing primer binding site




Sequencing with Index Reads




Single Index Reads

All Platforms



Dual Index Reads


MiSeq, HiSeq 2500, NovaSeq 6000



Dual indexed sequencing utilizes four sequencing reads

Dual Index Reads


iSeq 100, MiniSeq, NextSeq, HiSeq 3000/4000



Illumina Sequencing Workflow


1

Library Preparation




2

Cluster Generation



3

Sequencing









4

Data Analysis

```
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
5TAAGGCTAGGTTTCATGCTA
```

Analysis Overview

Analysis Type	Software	Outputs
Control Software		 Images, Intensities and Base Calls
Analysis Software		 Alignments, Variant Detection
Visualization Software		 Annotation, Filtering, Reports

Questions?

